

Antibodies Product Datasheet

Product Name:	AKAP7 Antibody
Catalog Number:	C48565
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AKAP7
Concentration:	1 mg/ml
Gene Name:	AKAP7 AKAP15 AKAP18
Protein Name:	A-kinase anchor protein 7 isoforms alpha and beta (AKAP-7 isoforms alpha and beta) (A-kinase anchor protein 18 kDa) (AKAP 18) (Protein kinase A-anchoring protein 7 isoforms alpha/beta) (PRKA7 isoforms alpha/beta)
Other Name:	
Human Swiss Prot No.:	O43687
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AKAP7 at Human
Storage and Stability:	Store at -20° C/1 year