

Antibodies Product Datasheet

Product Name:	EPHX2 Antibody
Catalog Number:	C48566
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHX2
Concentration:	1 mg/ml
Gene Name:	EPHX2
Protein Name:	Bifunctional epoxide hydrolase 2 [Includes: Cytosolic epoxide hydrolase 2 (CEH) (EC 3.3.2.10) (Epoxide hydratase) (Soluble epoxide hydrolase) (SEH)Lipid-phosphate phosphatase (EC 3.1.3.76)]
Other Name:	
Human Swiss Prot No.:	P34913
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EPHX2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year