

Antibodies Product Datasheet

Product Name:	ARL8A Antibody
Catalog Number:	C48567
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARL8A
Concentration:	1 mg/ml
Gene Name:	ARL8A ARL10B GIE2
Protein Name:	ADP-ribosylation factor-like protein 8A (ADP-ribosylation factor-like protein 10B) (Novel small G protein indispensable for equal chromosome segregation 2)
Other Name:	
Human Swiss Prot No.:	Q96BM9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARL8A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year