

Antibodies Product Datasheet

Product Name:	GCS1 Antibody
Catalog Number:	C48586
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GCS1
Concentration:	1 mg/ml
Gene Name:	MOGS GCS1
Protein Name:	Mannosyl-oligosaccharide glucosidase (EC 3.2.1.106) (Processing A-glucosidase I)
Other Name:	
Human Swiss Prot No.:	Q13724
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GCS1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year