

Antibodies Product Datasheet

Product Name:	KCC2/SLC12A5 Antibody
Catalog Number:	C48587
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KCC2/SLC12A5
Concentration:	1 mg/ml
Gene Name:	SLC12A5 KCC2 KIAA1176
Protein Name:	Solute carrier family 12 member 5 (Electroneutral potassium-chloride cotransporter 2) (K-Cl cotransporter 2) (hKCC2) (Neuronal K-Cl cotransporter)
Other Name:	
Human Swiss Prot No.:	Q9H2X9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KCC2/SLC12A5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year