

Antibodies Product Datasheet

Product Name:	HSPA13 Antibody
Catalog Number:	C48618
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HSPA13
Concentration:	1 mg/ml
Gene Name:	HSPA13 STCH
Protein Name:	Heat shock 70 kDa protein 13 (Microsomal stress-70 protein ATPase core) (Stress-70 protein chaperone microsome-associated 60 kDa protein)
Other Name:	
Human Swiss Prot No.:	P48723
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HSPA13 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year