

Antibodies Product Datasheet

Product Name:	DBNL Antibody
Catalog Number:	C48623
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DBNL
Concentration:	1 mg/ml
Gene Name:	DBNL CMAP SH3P7 PP5423
Protein Name:	Drebrin-like protein (Cervical SH3P7) (Cervical mucin-associated protein) (Drebrin-F) (HPK1-interacting protein of 55 kDa) (HIP-55) (SH3 domain-containing protein 7)
Other Name:	
Human Swiss Prot No.:	Q9UJU6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DBNL at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year