

Antibodies Product Datasheet

Product Name:	AADAT Antibody
Catalog Number:	C48631
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AADAT
Concentration:	1 mg/ml
Gene Name:	AADAT KAT2
Protein Name:	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial (KAT/AadAT) (2-aminoadipate aminotransferase) (2-aminoadipate transaminase) (EC 2.6.1.39) (Alpha-aminoadipate aminotransferase) (AadAT) (Kynurenine aminotransferase II) (Kynurenine--oxoglutarat
Other Name:	
Human Swiss Prot No.:	Q8N5Z0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AADAT at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year