

Antibodies Product Datasheet

Product Name:	KIAA0319L Antibody
Catalog Number:	C48641
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KIAA0319L
Concentration:	1 mg/ml
Gene Name:	KIAA0319L KIAA1837 PP791
Protein Name:	Dyslexia-associated protein KIAA0319-like protein
Other Name:	
Human Swiss Prot No.:	Q8IZA0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KIAA0319L at Human, Mouse
Storage and Stability:	Store at -20° C/1 year