

Antibodies Product Datasheet

Product Name:	DDHD2 Antibody
Catalog Number:	C48643
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DDHD2
Concentration:	1 mg/ml
Gene Name:	DDHD2 KIAA0725 SAMWD1
Protein Name:	Phospholipase DDHD2 (EC 3.1.1.-) (DDHD domain-containing protein 2) (SAM, WWE and DDHD domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	O94830
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DDHD2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year