

Antibodies Product Datasheet

Product Name:	ERCC6L Antibody
Catalog Number:	C48659
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ERCC6L
Concentration:	1 mg/ml
Gene Name:	ERCC6L PICH
Protein Name:	DNA excision repair protein ERCC-6-like (EC 3.6.4.12) (ATP-dependent helicase ERCC6-like) (PLK1-interacting checkpoint helicase) (Tumor antigen BJ-HCC-15)
Other Name:	
Human Swiss Prot No.:	Q2NKX8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ERCC6L at Human, Mouse
Storage and Stability:	Store at -20° C/1 year