

Antibodies Product Datasheet

Product Name:	BAG6 Antibody
Catalog Number:	C48677
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BAG6
Concentration:	1 mg/ml
Gene Name:	BAG6 BAT3 G3
Protein Name:	Large proline-rich protein BAG6 (BAG family molecular chaperone regulator 6) (BCL2-associated athanogene 6) (BAG-6) (BAG6) (HLA-B-associated transcript 3) (Protein G3) (Protein Scythe)
Other Name:	
Human Swiss Prot No.:	P46379
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BAG6 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year