

Antibodies Product Datasheet

Product Name:	DTX4 Antibody
Catalog Number:	C48681
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DTX4
Concentration:	1 mg/ml
Gene Name:	DTX4 KIAA0937 RNF155
Protein Name:	E3 ubiquitin-protein ligase DTX4 (EC 6.3.2.-) (Protein deltex-4) (Deltex4) (RING finger protein 155)
Other Name:	
Human Swiss Prot No.:	Q9Y2E6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DTX4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year