

Antibodies Product Datasheet

Product Name:	EIF3B Antibody
Catalog Number:	C48692
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF3B
Concentration:	1 mg/ml
Gene Name:	EIF3B EIF3S9
Protein Name:	Eukaryotic translation initiation factor 3 subunit B (eIF3b) (Eukaryotic translation initiation factor 3 subunit 9) (Prt1 homolog) (hPrt1) (eIF-3-eta) (eIF3 p110) (eIF3 p116)
Other Name:	
Human Swiss Prot No.:	P55884
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF3B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year