

Antibodies Product Datasheet

Product Name:	Kindlin 1 Antibody
Catalog Number:	C48700
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Kindlin 1
Concentration:	1 mg/ml
Gene Name:	FERMT1 C20orf42 KIND1 URP1
Protein Name:	Fermitin family homolog 1 (Kindlerin) (Kindlin syndrome protein) (Kindlin-1) (Unc-112-related protein 1)
Other Name:	
Human Swiss Prot No.:	Q9BQL6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Kindlin 1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year