

Antibodies Product Datasheet

Product Name:	GMAP-210 Antibody
Catalog Number:	C48715
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GMAP-210
Concentration:	1 mg/ml
Gene Name:	TRIP11 CEV14
Protein Name:	Thyroid receptor-interacting protein 11 (TR-interacting protein 11) (TRIP-11) (Clonal evolution-related gene on chromosome 14 protein) (Golgi-associated microtubule-binding protein 210) (GMAP-210) (Trip230)
Other Name:	
Human Swiss Prot No.:	Q15643
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GMAP-210 at Human
Storage and Stability:	Store at -20° C/1 year