

Antibodies Product Datasheet

Product Name:	LANCL1 Antibody
Catalog Number:	C48719
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LANCL1
Concentration:	1 mg/ml
Gene Name:	LANCL1 GPR69A
Protein Name:	LanC-like protein 1 (40 kDa erythrocyte membrane protein) (p40)
Other Name:	
Human Swiss Prot No.:	O43813
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LANCL1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year