

Antibodies Product Datasheet

Product Name:	EIF3G Antibody
Catalog Number:	C48723
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein. AA range:C-teral
Concentration:	1 mg/ml
Gene Name:	EIF3G EIF3S4
Protein Name:	Eukaryotic translation initiation factor 3 subunit G (eIF3g) (Eukaryotic translation initiation factor 3 RNA-binding subunit) (eIF-3 RNA-binding subunit) (Eukaryotic translation initiation factor 3 subunit 4) (eIF-3-delta) (eIF3 p42) (eIF3 p44)
Other Name:	
Human Swiss Prot No.:	O75821
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF3G at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year