

Antibodies Product Datasheet

Product Name:	APH1A Antibody
Catalog Number:	C48731
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human APH1A
Concentration:	1 mg/ml
Gene Name:	APH1A PSF CGI-78 UNQ579/PRO1141
Protein Name:	Gamma-secretase subunit APH-1A (APH-1a) (Aph-1alpha) (Presenilin-stabilization factor)
Other Name:	
Human Swiss Prot No.:	Q96BI3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of APH1A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year