

Antibodies Product Datasheet

Product Name:	IA-2/PTPRN Antibody
Catalog Number:	C48743
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IA-2/PTPRN
Concentration:	1 mg/ml
Gene Name:	PTPRN ICA3 ICA512
Protein Name:	Receptor-type tyrosine-protein phosphatase-like N (R-PTP-N) (Islet cell antigen 512) (ICA 512) (Islet cell autoantigen 3) (PTP IA-2)
Other Name:	
Human Swiss Prot No.:	Q16849
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IA-2/PTPRN at Human,Rat
Storage and Stability:	Store at -20° C/1 year