

Antibodies Product Datasheet

Product Name:	IRSp53 Antibody
Catalog Number:	C48746
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IRSp53
Concentration:	1 mg/ml
Gene Name:	BAIAP2
Protein Name:	Brain-specific angiogenesis inhibitor 1-associated protein 2 (BAI-associated protein 2) (BAI1-associated protein 2) (Protein BAP2) (Fas ligand-associated factor 3) (FLAF3) (Insulin receptor substrate p53/p58) (IRS-58) (IRSp53/58) (Insulin receptor substra
Other Name:	
Human Swiss Prot No.:	Q9UQB8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IRSp53 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year