

Antibodies Product Datasheet

Product Name:	HYAL2 Antibody
Catalog Number:	C48750
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HYAL2
Concentration:	1 mg/ml
Gene Name:	HYAL2 LUCA2
Protein Name:	Hyaluronidase-2 (Hyal-2) (EC 3.2.1.35) (Hyaluronoglucosaminidase-2) (Lung carcinoma protein 2) (LuCa-2)
Other Name:	
Human Swiss Prot No.:	Q12891
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HYAL2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year