

Antibodies Product Datasheet

Product Name:	ABCB6 Antibody
Catalog Number:	C48755
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ABCB6
Concentration:	1 mg/ml
Gene Name:	ABCB6 MTABC3 PRP UMAT
Protein Name:	ATP-binding cassette sub-family B member 6, mitochondrial (Mitochondrial ABC transporter 3) (Mt-ABC transporter 3) (P-glycoprotein-related protein) (Ubiquitously-expressed mammalian ABC half transporter)
Other Name:	
Human Swiss Prot No.:	Q9NP58
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ABCB6 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year