

Antibodies Product Datasheet

Product Name:	BDH2 Antibody
Catalog Number:	C48759
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BDH2
Concentration:	1 mg/ml
Gene Name:	BDH2 DHRS6 UNQ6308/PRO20933
Protein Name:	3-hydroxybutyrate dehydrogenase type 2 (EC 1.1.1.-) (EC 1.1.1.30) (Dehydrogenase/reductase SDR family member 6) (Oxidoreductase UCPA) (R-beta-hydroxybutyrate dehydrogenase)
Other Name:	
Human Swiss Prot No.:	Q9BUT1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BDH2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year