

Antibodies Product Datasheet

Product Name:	BPIL1 Antibody
Catalog Number:	C48767
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BPIL1
Concentration:	1 mg/ml
Gene Name:	BPIFB2 BPIL1 C20orf184 LPLUNC2 UNQ2489/PRO5776
Protein Name:	BPI fold-containing family B member 2 (Bactericidal/permeability-increasing protein-like 1) (BPI-like 1) (Long palate, lung and nasal epithelium carcinoma-associated protein 2) (RYSR)
Other Name:	
Human Swiss Prot No.:	Q8N4F0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BPIL1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year