

Antibodies Product Datasheet

Product Name:	KDM2A Antibody
Catalog Number:	C48777
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KDM2A
Concentration:	1 mg/ml
Gene Name:	KDM2A CXXC8 FBL7 FBXL11 JHDM1A KIAA1004
Protein Name:	Lysine-specific demethylase 2A (EC 1.14.11.27) (CXXC-type zinc finger protein 8) (F-box and leucine-rich repeat protein 11) (F-box protein FBL7) (F-box protein Lilina) (F-box/LRR-repeat protein 11) (JmjC domain-containing histone demethylation protein 1A)
Other Name:	
Human Swiss Prot No.:	Q9Y2K7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDM2A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year