

Antibodies Product Datasheet

Product Name:	HOP2 Antibody
Catalog Number:	C48786
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HOP2
Concentration:	1 mg/ml
Gene Name:	PSMC3IP HOP2 TBPIP
Protein Name:	Homologous-pairing protein 2 homolog (Nuclear receptor coactivator GT198) (PSMC3-interacting protein) (Proteasome 26S ATPase subunit 3-interacting protein) (Tat-binding protein 1-interacting protein) (TBP-1-interacting protein)
Other Name:	
Human Swiss Prot No.:	Q9P2W1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HOP2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year