

Antibodies Product Datasheet

Product Name:	DAB2IP Antibody
Catalog Number:	C48789
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DAB2IP
Concentration:	1 mg/ml
Gene Name:	DAB2IP AF9Q34 AIP1 KIAA1743
Protein Name:	Disabled homolog 2-interacting protein (DAB2 interaction protein) (DAB2-interacting protein) (ASK-interacting protein 1)
Other Name:	
Human Swiss Prot No.:	Q5VWQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DAB2IP at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year