

Antibodies Product Datasheet

Product Name:	ITPA Antibody
Catalog Number:	C48791
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ITPA
Concentration:	1 mg/ml
Gene Name:	ITPA C20orf37 My049 OK/SW-cl.9
Protein Name:	Inosine triphosphate pyrophosphatase (ITPase) (Inosine triphosphatase) (EC 3.6.1.19) (Non-canonical purine NTP pyrophosphatase) (Non-standard purine NTP pyrophosphatase) (Nucleoside-triphosphate diphosphatase) (Nucleoside-triphosphate pyrophosphatase) (NT
Other Name:	
Human Swiss Prot No.:	Q9BY32
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ITPA at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year