

Antibodies Product Datasheet

Product Name:	CTR9 Antibody
Catalog Number:	C48793
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CTR9.AA range:1150-1173
Concentration:	1 mg/ml
Gene Name:	CTR9 KIAA0155 SH2BP1
Protein Name:	RNA polymerase-associated protein CTR9 homolog (SH2 domain-binding protein 1)
Other Name:	
Human Swiss Prot No.:	Q6PD62
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CTR9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year