

Antibodies Product Datasheet

Product Name:	COPS8/COP9 Antibody
Catalog Number:	C48806
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COPS8/COP9
Concentration:	1 mg/ml
Gene Name:	COPS8 CSN8
Protein Name:	COP9 signalosome complex subunit 8 (SGN8) (Signalosome subunit 8) (COP9 homolog) (hCOP9) (JAB1-containing signalosome subunit 8)
Other Name:	
Human Swiss Prot No.:	Q99627
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COPS8/COP9 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year