

Antibodies Product Datasheet

Product Name:	HSD17B3 Antibody
Catalog Number:	C48815
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HSD17B3
Concentration:	1 mg/ml
Gene Name:	HSD17B3 EDH17B3
Protein Name:	Testosterone 17-beta-dehydrogenase 3 (EC 1.1.1.64) (17-beta-hydroxysteroid dehydrogenase type 3) (17-beta-HSD 3) (Testicular 17-beta-hydroxysteroid dehydrogenase)
Other Name:	
Human Swiss Prot No.:	P37058
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HSD17B3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year