

Antibodies Product Datasheet

Product Name:	LETM2 Antibody
Catalog Number:	C48816
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LETM2
Concentration:	1 mg/ml
Gene Name:	LETM2
Protein Name:	LETM1 domain-containing protein LETM2, mitochondrial (LETM1 and EF-hand domain-containing protein 2) (Leucine zipper-EF-hand-containing transmembrane protein 1-like)
Other Name:	
Human Swiss Prot No.:	Q2VYF4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LETM2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year