

Antibodies Product Datasheet

Product Name:	ATP11B Antibody
Catalog Number:	C48819
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ATP11B
Concentration:	1 mg/ml
Gene Name:	ATP11B ATP1F ATP1R KIAA0956
Protein Name:	Probable phospholipid-transporting ATPase IF (EC 3.6.3.1) (ATPase IR) (ATPase class VI type 11B)
Other Name:	
Human Swiss Prot No.:	Q9Y2G3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATP11B at Human
Storage and Stability:	Store at -20° C/1 year