

Antibodies Product Datasheet

Product Name:	ABP1 Antibody
Catalog Number:	C48821
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ABP1
Concentration:	1 mg/ml
Gene Name:	ABP1 AOC1 DAO1
Protein Name:	Amiloride-sensitive amine oxidase [copper-containing] (DAO) (Diamine oxidase) (EC 1.4.3.22) (Amiloride-binding protein) (ABP) (Histaminase) (Kidney amine oxidase) (KAO)
Other Name:	
Human Swiss Prot No.:	P19801
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ABP1 at Human,Rat
Storage and Stability:	Store at -20° C/1 year