

Antibodies Product Datasheet

Product Name:	ARIH2 Antibody
Catalog Number:	C48848
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARIH2
Concentration:	1 mg/ml
Gene Name:	ARIH2 ARI2 TRIAD1 HT005
Protein Name:	E3 ubiquitin-protein ligase ARIH2 (ARI-2) (Protein ariadne-2 homolog) (EC 6.3.2.-) (Triad1 protein)
Other Name:	
Human Swiss Prot No.:	O95376
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARIH2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year