

Antibodies Product Datasheet

Product Name:	Adenylosuccinate lyase Antibody
Catalog Number:	C48857
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Adenylosuccinate lyase
Concentration:	1 mg/ml
Gene Name:	ADSL AMPS
Protein Name:	Adenylosuccinate lyase (ASL) (EC 4.3.2.2) (Adenylosuccinase) (ASase)
Other Name:	
Human Swiss Prot No.:	P30566
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Adenylosuccinate lyase at Human, Mouse
Storage and Stability:	Store at -20° C/1 year