

Antibodies Product Datasheet

Product Name:	EXOSC4 Antibody
Catalog Number:	C48864
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EXOSC4
Concentration:	1 mg/ml
Gene Name:	EXOSC4 RRP41 SKI6
Protein Name:	Exosome complex component RRP41 (Exosome component 4) (Ribosomal RNA-processing protein 41) (p12A)
Other Name:	
Human Swiss Prot No.:	Q9NPD3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EXOSC4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year