

Antibodies Product Datasheet

Product Name:	AGL Antibody
Catalog Number:	C48865
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AGL
Concentration:	1 mg/ml
Gene Name:	AGL GDE
Protein Name:	Glycogen debranching enzyme (Glycogen debrancher) [Includes: 4-alpha-glucanotransferase (EC 2.4.1.25) (Oligo-1,4-1,4-glucantransferase)Amylo-alpha-1,6-glucosidase (Amylo-1,6-glucosidase) (EC 3.2.1.33) (Dextrin 6-alpha-D-glucosidase)]
Other Name:	
Human Swiss Prot No.:	P35573
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AGL at Human
Storage and Stability:	Store at -20° C/1 year