

Antibodies Product Datasheet

Product Name:	ASIC2 Antibody
Catalog Number:	C48901
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASIC2
Concentration:	1 mg/ml
Gene Name:	ASIC2 ACCN ACCN1 BNAC1 MDEG
Protein Name:	Acid-sensing ion channel 2 (ASIC2) (Amiloride-sensitive brain sodium channel) (Amiloride-sensitive cation channel 1, neuronal) (Amiloride-sensitive cation channel neuronal 1) (Brain sodium channel 1) (BNC1) (BNAC1) (Mammalian degenerin homolog)
Other Name:	
Human Swiss Prot No.:	Q16515
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASIC2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year