

Antibodies Product Datasheet

Product Name:	ACOT7 Antibody
Catalog Number:	C48903
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ACOT7
Concentration:	1 mg/ml
Gene Name:	ACOT7 BACH
Protein Name:	Cytosolic acyl coenzyme A thioester hydrolase (EC 3.1.2.2) (Acyl-CoA thioesterase 7) (Brain acyl-CoA hydrolase) (BACH) (CTE-IIa) (CTE-II) (Long chain acyl-CoA thioester hydrolase)
Other Name:	
Human Swiss Prot No.:	O00154
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ACOT7 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year