

Antibodies Product Datasheet

Product Name:	GIYD2 Antibody
Catalog Number:	C48908
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GIYD2
Concentration:	1 mg/ml
Gene Name:	SLX1A GIYD1 SLX1; SLX1B GIYD2 SLX1
Protein Name:	Structure-specific endonuclease subunit SLX1 (EC 3.6.1.-) (GIY-YIG domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q9BQ83
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GIYD2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year