

Antibodies Product Datasheet

Product Name:	CLCA2 Antibody
Catalog Number:	C48935
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CLCA2
Concentration:	1 mg/ml
Gene Name:	CLCA2 CACC3
Protein Name:	Calcium-activated chloride channel regulator 2 (Calcium-activated chloride channel family member 2) (hCLCA2) (Calcium-activated chloride channel protein 3) (CaCC-3) (hCaCC-3) [Cleaved into: Calcium-activated chloride channel regulator 2, 109 kDa formCalci
Other Name:	
Human Swiss Prot No.:	Q9UQC9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CLCA2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year