

Antibodies Product Datasheet

Product Name:	FUT7 Antibody
Catalog Number:	C48942
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FUT7
Concentration:	1 mg/ml
Gene Name:	FUT7
Protein Name:	Alpha-(1,3)-fucosyltransferase (EC 2.4.1.-) (Fucosyltransferase 7) (Fucosyltransferase VII) (Fuc-TVII) (FucT-VII) (Galactoside 3-L-fucosyltransferase) (Selectin ligand synthase)
Other Name:	
Human Swiss Prot No.:	Q11130
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FUT7 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year