

Antibodies Product Datasheet

Product Name:	BANP Antibody
Catalog Number:	C48950
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BANP
Concentration:	1 mg/ml
Gene Name:	BANP BEND1 SMAR1
Protein Name:	Protein BANP (BEN domain-containing protein 1) (Btg3-associated nuclear protein) (Scaffold/matrix-associated region-1-binding protein)
Other Name:	
Human Swiss Prot No.:	Q8N9N5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BANP at Human, Mouse
Storage and Stability:	Store at -20° C/1 year