

Antibodies Product Datasheet

Product Name:	FAM119A Antibody
Catalog Number:	C48960
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FAM119A
Concentration:	1 mg/ml
Gene Name:	METTL21A FAM119A HCA557B
Protein Name:	Methyltransferase-like protein 21A (EC 2.1.1.-) (Hepatocellular carcinoma-associated antigen 557b)
Other Name:	
Human Swiss Prot No.:	Q8WXB1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FAM119A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year