

Antibodies Product Datasheet

Product Name:	Galectin 10 Antibody
Catalog Number:	C48963
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Galectin 10
Concentration:	1 mg/ml
Gene Name:	CLC
Protein Name:	Eosinophil lysophospholipase (EC 3.1.1.5) (Charcot-Leyden crystal protein) (CLC) (Galectin-10) (Gal-10) (Lysolecithin acylhydrolase)
Other Name:	
Human Swiss Prot No.:	Q05315
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Galectin 10 at Human
Storage and Stability:	Store at -20° C/1 year