

Antibodies Product Datasheet

Product Name:	LSM4 Antibody
Catalog Number:	C48965
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LSM4
Concentration:	1 mg/ml
Gene Name:	LSM4
Protein Name:	U6 snRNA-associated Sm-like protein LSm4 (Glycine-rich protein) (GRP)
Other Name:	
Human Swiss Prot No.:	Q9Y4Z0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LSM4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year