

Antibodies Product Datasheet

Product Name:	CAPRIN2 Antibody
Catalog Number:	C48972
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CAPRIN2
Concentration:	1 mg/ml
Gene Name:	CAPRIN2 C1QDC1 EEG1 KIAA1873 RNG140
Protein Name:	Caprin-2 (C1q domain-containing protein 1) (Cytoplasmic activation/proliferation-associated protein 2) (Gastric cancer multidrug resistance-associated protein) (Protein EEG-1) (RNA granule protein 140)
Other Name:	
Human Swiss Prot No.:	Q6IMN6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CAPRIN2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year