

Antibodies Product Datasheet

Product Name:	BRD8 Antibody
Catalog Number:	C48973
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BRD8
Concentration:	1 mg/ml
Gene Name:	BRD8 SMAP SMAP2
Protein Name:	Bromodomain-containing protein 8 (Skeletal muscle abundant protein) (Skeletal muscle abundant protein 2) (Thyroid hormone receptor coactivating protein of 120 kDa) (TrCP120) (p120)
Other Name:	
Human Swiss Prot No.:	Q9H0E9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRD8 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year